CHAMOT

**Annexin V-iFluor 488/PI Apoptosis Detection Reagent 细胞凋亡
检测试剂盒**

CM001-50D


CM001-100D

CHAMOT

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Specialize In Cytokines

CONTENT

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Annexin V-iFluor 488/PI Apoptosis Detection Reagent 细胞凋亡检测试剂盒

产品编号	规格
CM001-50D	50 rxn
CM001-100D	100 rxn

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- Effect of Metformin Nanoparticle-Mediated Thioredoxin Interacting Protein Expression on Oxaliplatin-Induced Peripheral Neuropathy. *J Nanosci Nanotechnol.* 2020 Oct 1;20(10):6123-6132. (IF 1.42)
- Platelets are highly efficient and efficacious carriers for tumor-targeted nano-drug delivery. *Drug Delivery.* 2022, VOL. 29, NO. 1, 937–949. (IF 6.419)
- A nanoreactor boosts chemodynamic therapy and ferroptosis for synergistic cancer therapy using molecular amplifier dihydroartemisinin. *J Nanobiotechnology.* 2022 May 14;20(1):230. (IF 9.464)
- Dihydroartemisinin remodels macrophage into an M1 phenotype via ferroptosis-mediated DNA damage. *Front Pharmacol.* 2022 Aug 11;13:949835. (IF 5.81)
- Targeted photodynamic therapy of glioblastoma mediated by platelets with photo-controlled release property. *Biomaterials.* 2022 Sep 30;290:121833. (IF 15.304)
- Antagonizing apolipoprotein J chaperone promotes proteasomal degradation of mTOR and relieves hepatic lipid deposition. *Hepatology.* 2023 Jan 3. (IF 17.298)
- Chlorin e6-induced photodynamic effect facilitates immunogenic cell death of lung cancer as a result of oxidative endoplasmic reticulum stress and DNA damage. *Int Immunopharmacol.* 2023 Feb;115:109661. (IF 5.714)
- Dihydroartemisinin elicits immunogenic death through ferroptosis-triggered ER stress and DNA damage for lung cancer immunotherapy. *Phytomedicine.* 2023 Apr;112:154682. (5.34)
- Laser-triggered intelligent drug delivery and anti-cancer photodynamic therapy using platelets as the vehicle. *Platelets.* 2023 Dec;34(1):2166677. (IF 4.236)
- Up-regulation of ABCG2 by MYBL2 deletion drives Chlorin e6-mediated photodynamic therapy resistance in colorectal cancer. *Photodiagnosis Photodyn Ther.* 2023 Apr 6;103558. (IF 3.577)
- Sonodynamic therapy of glioblastoma mediated by platelets with ultrasound-triggered drug release. *Drug Deliv.* 2023 Dec;30(1):2219429. (IF 6.819)
- MiR-934 Exacerbates Malignancy of Gastric Cancer Cells by Targeting ZFP36. *Iran J Public Health.* 2023;52(8):1720-1729. (IF 1.4)
- A Dihydroartemisinin-Loaded Nanoreactor Motivates Anti-Cancer Immunotherapy by Synergy-Induced Ferroptosis to Activate Cgas/STING for Reprogramming of Macrophage. *Adv Healthc Mater.* 2023 Aug 11;e2301561. (IF 11.092)
- Anti-lung cancer synergy of low-dose doxorubicin and PD-L1 blocker co-delivered via mild photothermia-responsive black phosphorus. *Drug Deliv Transl Res.* 2024 Apr 10. (IF 5.4)
- Dihydroartemisinin-driven TOM70 inhibition leads to mitochondrial destabilization to induce pyroptosis against lung cancer. *Phytother Res.* 2024 May 18. (IF 6.1)
- Targeted Sonodynamic Therapy Induces Tumor Cell Quasi-Immunogenic Ferroptosis and Macrophage Immunostimulatory Autophagy in

Glioblastoma. Biomaterials. Volume 315, April 2025, 122913. (IF 12.8)

- Suppression of NRF2 by photodynamic action to enhance bortezomib-triggered DNA damage for synergistic colorectal cancer therapy. J

Photochem Photobiol B. 2025 Jun;267:113162. (IF 4.5)

- Accelerated differentiation of photothermal effect-induced mesenchymal stem cells regulated by activating HSP90-autophagy boosts wound repair.

Int J Pharm. 2025 May 15;676:125631. (IF 5.77)

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491