

CHAMOT

Mouse anti-GST tag mAb

CM005-0.1A1

CM005-0.5A1

CM005-1A1

CM005-5A1

CM005-10A1



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Mouse anti-GST tag mAb

产品编号	规格
CM005-0.1A1	100 µg
CM005-0.5A1	500 µg
CM005-1A1	1 mg
CM005-5A1	5 mg
CM005-10A1	10 mg

产品简介

产品形式	Liquid
产品应用	ELISA,WB,IP,IFA,FACS
产品描述	Mouse anti-GST tag mAb is specific for native as well as denatured-reduced forms of GST (glutathione S-transferase) tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and FACS.
反应性	Recognizes glutathione S-transferase tag on either the N- or C-terminus of recombinant fusion protein.
克隆性	Monoclonal, clone 89-3
免疫原	Recombinant glutathione S-transferase
浓度	1 mg/mL
标记	Unconjugated
产品组成	Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.
亚型	IgG2a

产品使用

Recommended dilution factor:

ELISA: 1:5000-20000

WB:1:1000-5000

IP: 1:500-2000

IFA: 1:500-1000

FACS: 1:500-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only
Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com
Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491