


# CHAMOT

COVID-19 Human Control Antibody Set (NP)(IgG/IgA/IgM)

CM013-0.02A1CoV



CHAMOT

喬默生物技術(上海)有限公司  
CHAMOT BIOTECHNOLOGY CO., LTD.

# CONTENT

1 产品简介

2 产品储存

3 产品组成

4 产品应用

## COVID-19 Human Control Antibody Set (NP)(IgG/IgA/IgM)

编号:	CM013-0.02A1CoV	规格:	20 µg
类别:	一抗	应用:	ELISA /WB/ IFA/FACS

### 产品简介

产品描述:	Human anti-SARS-CoV & CoV-2 NP antibody recognizes human SARS Coronavirus Nucleocapsid protein (NP). Coronavirus NP localizes to the cytoplasm and the nucleolus in both Virus-like particle (VLP) infected primary cells and in cells transfected with NP plasmid. NP has abundant expression in coronavirus, and is a highly immunogenic phosphoprotein. NP is also conserved in sequence. Due to these characteristics above, the NP is an ideal marker for diagnosis.
类型:	重组单克隆抗体
浓度:	取决于批次号
免疫原:	N/A
宿主:	Human
反应种属:	Virus
偶联标记:	N/A
物理状态:	液体

### 产品储存

储存:	4°C保存六个月; -20°C可长期保存, 避免反复冻融。
运输:	蓝冰。

### 产品组成

含 pH7.4 PBS 溶液。

WB=Western Blot; IP=Immunoprecipitation; IF=Immunofluorescence; IHC=Immunohistochemistry;  
FACS=Fluorescence activated Cell Sorting; FA=Functional Assay