

CHAMOT

Recombinant IL-10, Human, GMP

CM012-100HPG
CM012-1000HPG



CHAMOT

喬默生物技術(上海)有限公司
CHAMOT BIOTECHNOLOGY CO., LTD.

CONTENT

1 产品简介

2 产品组成

3 产品储存

4 产品使用


5 产品规范

Recombinant IL-10 (Interleukin-10), Human, GMP

编号:	CM012-100HPG CM012-1000HPG	规格:	100 µg 1 mg
类别:	重组蛋白	应用:	Functional Assay

产品简介

描述: Interleukin-10 (IL-10), also known as human cytokine synthesis inhibitory factor (CSIF), is an anti-inflammatory cytokine. In humans, IL-10 is encoded by the IL-10 gene. IL-10 signals through a receptor complex consisting of two IL-10 receptor-1 and two IL-10 receptor-2 proteins. Consequently, the functional receptor consists of four IL-10 receptor molecules. IL-10 binding induces STAT3 signaling via the phosphorylation of the cytoplasmic tails of IL-10 receptor-1 + IL-10 receptor-2 by JAK1 and Tyk2 respectively.



来源: *Escherichia coli*.
Animal-free reagent and laboratory.
Manufactured and tested under GMP guideline.

SDS-PAGE analysis of recombinant IL-10, Human, GMP

纯度: >98% as determined by SDS-PAGE analysis.

生物学活性: Measure by its ability to induce MC/9-2 cells proliferation. The ED₅₀ for this effect is <1 ng/mL. The specific activity of recombinant human IL-10 is approximately >1 x10⁶ IU/ mg.

内毒素检测: <0.05 EU per 1 µg of the protein by the LAL method.

蛋白序列: MSPGQGTOSENSCTHFPGNLPNMLRDLRDA
FSRVKTFQMKDQLDNLLKELLEDFKGYLGC
QALSEMIQFYLEEVM POAENQDPDIKAHVNS
LGENLKTLLR RLRRCHRFLPCENKSKAVEQVK
NAFNKLQEKGIYKAMSEFDIFINYIEAYMTMKI



RN with polyhistidine tag at the C - terminus

产品组成

成分: 从含有 1X PBS, pH 8.0 溶液中冻干的蛋白质, 无任何动物源性成分.

产品储存

储存: -20°C冻干粉保存 1 年.
运输: 蓝冰.

产品使用

- 1、开盖前, 建议3000-3500rpm离心5min。
- 2、推荐使用无菌水重悬冻干粉, 溶液浓度不低于500µg/mL, 并室温静置至少20min以充分溶解。勿涡旋剧烈振荡。
- 3、重悬后的溶液, 2-8无菌保存不超过1周。
- 4、如需长期保存, 推荐使用无菌的含载体蛋白(如0.1%BSA、10%FBS 或5%HSA)的溶液进一步稀释(不低于10ug/mL)后分装保存, -20至-80无菌保存3到6个月。无血清实验需求时, 可更换为 5%海藻糖溶液作为载体。避免反复冻融。

产品规范

CHAMOT GMP recombinant proteins are manufactured in ISO 13485:2016 and GMP certified facility. The processes include:

- ★Testing and traceability of raw material
- ★Records of the maintenance and equipment calibration
- ★Personnel training records
- ★Batch-to-batch consistency
- ★Documentation of QA control and process changes
- ★Manufactured and tested under an ISO 13485:2016 certified quality management system
- ★Stability monitor of product shelf-life

WB=Western Blot; IP=Immunoprecipitation; IF=Immunofluorescence; IHC=Immunohistochemistry;
FACS=Fluorescence activated Cell Sorting; FA=Functional Assay

CHAMOT Biotechnology (Shanghai) Co., Ltd.

For Research Use Only. Not for use in diagnostic or therapeutic procedures