

CHAMOT

Goat anti-Mouse IgG (H+L)

CM002-0.1A2

CM002-0.5A2

CM002-1A2

CM002-5A2

CM002-10A2



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品储存/运输

Goat anti-Mouse IgG (H+L)

产品编号	规格
CM002-0.1A2	100 µg
CM002-0.5A2	500 µg
CM002-1A2	1 mg
CM002-5A2	5 mg
CM002-10A2	10 mg

产品简介

产品形式	Liquid
产品描述	Goat Anti-Mouse IgG (H+L) Antibody is purified from goat, specifically react with mouse IgG, heavy and light chain.
反应性	Mouse (Minor cross-reactivity to rat)
克隆性	Polyclonal
免疫原	Purified mouse IgG
浓度	5 mg/mL
标记	Unconjugated
产品组成	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
亚型	IgG

产品储存/运输

储存	This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
运输	The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491