

CHAMOT

Goat anti-Mouse IgG-iFluor 405

CM011-0.1A2F

CM011-0.5A2F

CM011-1A2F

CM011-5A2F

CM011-10A2F



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

Goat anti-Mouse IgG-iFluor 405

产品编号	规格
CM011-0.1A2F	100 μL
CM011-0.5A2F	500 μL
CM011-1A2F	1 mL
CM011-5A2F	5 mL
CM011-10A2F	10 mL

产品简介

产品形式	Liquid
产品应用	IF,FACS
产品描述	iFluor 405 is a blue-fluorescent dye that can be excited by the 405 nm laser. LeadgeneGoat Anti-Mouse IgG (H+L) Antibody (iFluor 405) is conjugated with high-activity iFluor 405.
反应性	Mouse (Minor cross-reactivity to rat)
克隆性	Polyclonal
免疫原	Purified mouse IgG
浓度	0.5 mg/mL
标记	iFluor 405
产品组成	Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.
亚型	IgG

产品使用

Recommended dilution factor:

IF: 1:200-1000

FACS:1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存

This product is stable after storage at: 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.

运输

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491