

CHAMOT

Elisa试剂盒配套缓冲液 Matched Diluent

N/A-1

N/A-2

N/A-3



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品储存/运输

Elisa试剂盒配套缓冲液Matched Diluent

产品编号	规格
N/A-1	5 Plates
N/A-2	10 Plates
N/A-3	15 Plates

产品简介

别称	High Mobility Group Protein 1,HMG1
产品形式	Liquid
产品应用	ELISA
产品组成	包被液、样本/标准品稀释液、检测抗体稀释液

产品储存/运输

储存	4°C可长期保存
运输	Blue Ice

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com

Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491