

# CHAMOT

流感病毒A抗体 / Anti-Influenza A virus NP Antibody [clone IA03]


CM006-100AI



CHAMOT

乔默®生物

Specialize In Cytokines



# CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

## 流感病毒A抗体 / Anti-Influenza A virus NP Antibody [clone IA03]

产品编号	规格
CM006-100AI	100 µg

### 产品简介

产品形式	Liquid
产品应用	ELISA
产品描述	Anti-Influenza A virus NP Antibody only recognizes nucleocapsid protein (NP) of influenza A viruses but not influenza B viruses. Influenza A viruses are RNA viruses, and their subtypes are labeled according to an H number (for the type of hemagglutinin) and an N number (for the type of neuraminidase).
反应性	Influenza A virus Nucleocapsid protein
克隆性	Monoclonal, clone IA03
免疫原	Recombinant Influenza A virus nucleocapsid protein
浓度	1 mg/mL
标记	Unconjugated
产品组成	Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.
亚型	IgG2a

### 产品使用

Recognize Influenza A virus NP in Lateral Flow, when clone IA03 antibody was paired with Anti-Influenza A virus NP Antibody [clone IA01].

Recommended dilution factor:

ELISA: 1:5000-20000

### 产品储存/运输

#### 储存

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

## 运输

This product is stable after storage at: 2-8°C for 2 weeks under sterile conditions from date of receipt. -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

For Research Use or Further Manufacturing Only

Chamot Biotechnology(Shanghai) Co., Ltd. [www.chamot-bio.com](http://www.chamot-bio.com)

Tel: 021-51880030 Mail: [info@chamot-bio.com](mailto:info@chamot-bio.com) QQ: 864920491