

CHAMOT

呼吸道合胞病毒抗体 / Anti-Respiratory Syncytial Virus Antibody


CM001-100ARV



CHAMOT

乔默®生物

Specialize In Cytokines



CONTENT

1 产品简介

2 产品使用

3 产品储存/运输

4 实验数据展示

呼吸道合胞病毒抗体 / Anti-Respiratory Syncytial Virus Antibody

产品编号	规格
CM001-100ARV	100 ug

产品简介

别称	RSV
产品形式	Liquid
产品应用	ELISA,WB,IFA
反应性	recognizes the fusion protein in pre?and post fusion conformation. Reacts with RSV A2 , RSV B 9320, RSV Long.
克隆性	Monoclonal
免疫原	Human respiratory syncytial virus
浓度	1 mg/mL
标记	Unconjugated
产品组成	1×Phosphate Buffered Saline, pH 7.4.

产品使用

Recommended dilution factor:

ELISA: 1:5000-10000



WB: 1:5000-10000

IFA: 1:5000-10000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

产品储存/运输

储存	This product is stable after storage at:2-8°C for 6 months from date of receipt.-80°C for 24 months from date of receipt.Avoid repeated freeze/thaw cycles.Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
运输	Blue Ice



For Research Use or Further Manufacturing Only
Chamot Biotechnology(Shanghai) Co., Ltd. www.chamot-bio.com
Tel: 021-51880030 Mail: info@chamot-bio.com QQ: 864920491