

CHAMOT

Recombinant Human IL-18BP

CM194-5HP

CM194-20HP

CM194-100HP

CM194-500HP

CM194-1000HP



CHAMOT

喬默生物技術(上海)有限公司
CHAMOT BIOTECHNOLOGY CO., LTD.

CONTENT

1 产品简介

2 产品组成

3 产品储存


4 产品使用

Recombinant Human IL-18BP

编号:	CM194-5HP CM194-20HP CM194-100HP CM194-500HP CM194-1000HP	规格:	5 µg 20 µg 100 µg 500 µg 1 mg
类别:	重组蛋白	应用:	Functional Assay

产品简介

描述:	IL18BP, specifically binds to IL-18, suppressing its activity. IL18BP lacks a transmembrane domain and is prominently expressed in the heart, lung, placenta, and spleen. Additionally, IL18BP regulates the early TH1 cytokine response by reducing IFN-gamma production, moderating T-helper type 1 immune responses.
来源:	HEK293
纯度:	>95% as determined by SDS-PAGE.
生物学活性:	Testing in process
内毒素检测:	< 0.1 EU per 1 µg of the protein by the LAL method.
蛋白序列:	A DNA sequence encoding Human IL-18BP Protein (#O95998)(Met1-Gly194) was expressed with polyhistidine tag at the C-terminus.



SDS- PAGE analysis of recombinant human IL-18BP

产品组成

成分：从含有 1X PBS, pH 7.4 溶液中冻干的蛋白质。

产品储存/运输

产品形式	储存温度	储存时间
冻干粉	-20℃至-80℃	自收到之日起1年
重悬液 (初始)	2℃至8℃	不超过1周
重悬液 (经稀释)	-20℃至-80℃	3到6个月

运输方式：蓝冰

产品使用

- 1、开盖前，建议3000-3500rpm离心5 min。
- 2、推荐使用无菌水重悬冻干粉，溶液浓度不低于100μg/mL，不高于1mg/mL，并室温静置至少20min以充分溶解。勿涡旋剧烈振荡。
- 3、重悬后的溶液，2-8℃无菌保存不超过1周。
- 4、如需长期保存，推荐使用无菌的含载体蛋白（如0.1 %BSA、10 %FBS 或5 %HSA）的溶液进一步稀释（不低于10ug/mL）后分装保存，-20℃至 -80℃无菌保存3到6个月。无血清实验需求时，可更换为 5%海藻糖溶液作为载体。避免反复冻融。

WB= Western Blot; IP= Immunoprecipitation; IF= Immunofluorescence; IHC= Immunohistochemistry;
FACS= Fluorescence activated Cell Sorting; FA= Functional Assay