

CHAMOT

Recombinant Mouse BMP-4

CM049-5MP

CM049-20MP

CM049-100MP

CM049-500MP

CM049-1000MP



CHAMOT

喬默生物技術(上海)有限公司
CHAMOT BIOTECHNOLOGY CO., LTD.

CONTENT

1 产品简介

2 产品组成

3 产品储存

4 产品使用

Recombinant Mouse BMP-4 (Bone morphogenetic protein-4)

编号:	CM049-5MP CM049-20MP CM049-100MP CM049-500MP CM049-1000MP	规格:	5 µg 20 µg 100 µg 500 µg 1 mg
类别:	重组蛋白	应用:	Functional Assay

产品简介

描述: Bone morphogenetic protein 4 is a protein that in humans is encoded by BMP4 gene. BMP4 is found on chromosome 14q22-q23. BMP4 is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The superfamily includes large families of growth and differentiation factors. BMP4 is highly conserved evolutionarily. BMP4 is found in early embryonic development in the ventral marginal zone and in the eye, heart blood and otic vesicle.


来源: *Escherichia coli*

纯度: >98% as determined by SDS-PAGE.

生物学活性: Measure by its ability to induce alkaline phosphatase production by ATDC5 cells. The ED₅₀ for this effect is < 10 ng/mL. The specific activity of recombinant mouse BMP-4 is > 1 x 10⁵ IU/mg.

内毒素检测: < 0.1 EU per 1 µg of the protein by the LAL method.

蛋白序列: MKKNKNCRRHSLYVDFSDVGWND
WIVAPPGYQAFYCHGDCPFPLADHL
NSTNHAIVQTLVNSVNSSIPKACCV
TELSAISMLYLDEYDKVVLKNYQEMV
VEGCGCR with polyhistidine tag at the C-terminus



SDS-PAGE analysis of recombinant mouse BMP-4

产品组成

成分： 从含有 1X PBS, pH 8.0溶液中冻干的蛋白质.

产品储存/运输

产品形式	储存温度	储存时间
冻干粉	-20℃至-80℃	自收到之日起1年
重悬液（初始）	2℃至8℃	不超过1周
重悬液（经稀释）	-20℃至-80℃	3到6个月

运输方式： 蓝冰

产品使用

- 1、开盖前，建议3000-3500rpm离心5min。
- 2、推荐使用无菌水重悬冻干粉，溶液浓度不低于100μg/mL，不高于1mg/mL，并室温静置至少20min以充分溶解。勿涡旋剧烈振荡。
- 3、重悬后的溶液，2-8℃无菌保存不超过1周。
- 4、如需长期保存，推荐使用无菌的含载体蛋白（如0.1%BSA、10%FBS 或5%HSA）的溶液进一步稀释（不低于10ug/mL）后分装保存，-20℃至-80℃无菌保存3到6个月。无血清实验需求时，可更换为 5%海藻糖溶液作为载体。避免反复冻融。

WB= Western Blot; IP= Immunoprecipitation; IF= Immunofluorescence; IHC= Immunohistochemistry
 ;FACS= Fluorescence activated Cell Sorting; FA= Functional Assay