

CHAMOT

AccuBlot™ anti-Mouse IgG-HRP

CM009-0.1A2H

CM009-0.5A2H

CM009-1A2H



CHAMOT

喬默生物技術(上海)有限公司
CHAMOT BIOTECHNOLOGY CO., LTD.

CONTENT

1 产品简介

2 产品储存

3 产品组成

4 产品应用



AccuBlot™ anti-Mouse IgG-HRP

编号	CM009-0.1A2H CM009-0.5A2H CM009-1A2H	规格	100 µL 500 µL 1 mL
类别	酶标二抗	应用	WB/ELISA

产品简介

产品描述： AccuBlot™ Anti-Mouse IgG Antibody (HRP) is conjugated with high-activity HRP. In coimmunoprecipitation (co-IP) assay, denatured IgG from IP samples may show interference and result in immunoblotting at ~55 kDa (heavy chain) and ~25 kDa (light chain) by using traditional secondary antibodies. The products of AccuBlot™ series are suitable for native IgG detection.

来源：大鼠
反应种属：小鼠
抗体类型：单克隆抗体
免疫原：小鼠 IgG
特异性：IgG2b
标记：HRP
浓度：Lot-dependent

产品储存

储存：4℃ 两周
-20℃/-80℃ 1年
避免反复冻融
运输：蓝冰

产品组成

PBS 溶液 · 含 1%BSA · 0.03% ProClin 300和50% glycerol · pH 7.4

产品应用

推荐稀释比：

WB 1:2500-1:10000
ELISA 1:1000-1:3000

WB=Western Blot; IP=Immunoprecipitation; IF=Immunofluorescence; IHC=Immunohistochemistry;
FACS=Fluorescence activated Cell Sorting; FA=Functional Assay