

CHAMOT

Mouse anti-DYKDDDDDK mAb

CM006-0.1A1

CM006-0.5A1

CM006-1A1



CHAMOT

喬默生物技術(上海)有限公司
CHAMOT BIOTECHNOLOGY CO., LTD.

CONTENT

1 产品简介

2 产品储存

3 产品组成

4 产品应用

Mouse anti-DYKDDDDK mAb

编号	CM006-0.1A1 CM006-0.5A1 CM006-1A1	规格	100 µg/vial 500 µg/vial 1 mg/vial
类别	标签抗体	应用	ELISA/ WB/ IFA/ IHC/ IP

产品简介

描述：	Mouse anti-DYKDDDDK tag mAb is specific to DYKDDDDK(FLAG) tag which can be applied to several immunoassays including ELISA, WB, IP, IFA and IHC. This antibody has been identified in several sources (bacteria, yeast and mammalian cells) of FLAG tag fusion protein.
来源：	小鼠
类型：	单抗, clone108-21
同种型：	IgG2b kappa
反应性：	识别DYKDDDDK蛋白序列 (N端或C端)
免疫原：	多肽DYKDDDDK (FLAG)
浓度：	0.5 mg/mL
纯度：	>95%(SDS-PAGE)

产品储存

储存：	-20 °C 可长期保存 避免反复冻融
运输：	蓝冰

产品组成

组成：	PBS溶液, 含50% glycerol 和 0.02% sodium azide
-----	---

产品应用

推荐稀释比及终浓度：

Application	Dilution factor	Final concentration
ELISA	1:5000-20000	0.025-0.1 µg/mL
WB	1:1000-10000	0.05-0.5 µg/mL
IFA	1:500-2000	0.25-1 µg/mL
IHC	1:500-2000	0.25-1 µg/mL
IP	1:1000	0.5 µg/mL

WB=Western Blot; IP=Immunoprecipitation; IF=Immunofluorescence; IHC=Immunohistochemistry;
FACS=Fluorescence activated Cell Sorting; FA=Functional Assay